

Revision date: 2/25/2022

N109 W13300 ELLSWORTH DRIVE GERMANTOWN, WI 53022 262-253-5900 FAX 262-253-5919

DESCRIPTION:

Resinlab[®] *UR3010RC Clear* is a REACH compliant reformulation of UR3010 Clear.

UR3010RC Clear is a two-component, room temperature curing polyurethane potting system. It is designed to cure completely at room temperature. *UR3010RC Clear* provides excellent environmental protection giving very good resistance to water, salt spray, inorganic acids and bases, and most organic solvents, while maintaining adhesion to various metals, plastics and other common assembly materials.

This product was formulated to a convenient mixing ratio of 1 to 1 by volume. It is well suited to automated meter mix operations. *UR3010RC Clear* is a DOT non-hazardous material and contains no TDI, MOCA, or other reportable substances. Part B is moisture sensitive and must be kept away from atmospheric moisture during storage. After opening a sealed container, thoroughly purge the remaining air space with dry nitrogen (or equivalent) before closing the container.

TYPICAL PROPERTIES:

All properties given are at 25 °C unless otherwise noted.

All properties given are at 25°C unless otherwise noted.		
Property:	Value:	Test Method or Source:
Color	Clear	Visual
Mix Ratio	Part A to Part B	Calculated
Mix Ratio by weight	0.91 to 1	
Mix Ratio by volume	1 to 1	
Cure Schedule	24-48 hrs @ 25 °C	
	2 hrs @ 65 °C	
	30 min @ 100 °C	
Viscosity - Part A	750 cP	TA HR20 Rheometer 25mm parallel plate @
Viscosity - Part B	230 cP	1/s DCV6100723
Viscosity - Mixed	500 cP	
Specific Gravity - Part A	0.96	Calculated
Specific Gravity - Part B	1.07	
Specific Gravity - Mixed	1.02	
Pot Life defined as the time it takes for	10 minutes	TA HR20 Rheometer parallel plate 25mm @
initial mixed viscosity to double		1/s DCV6100723
Gel Time 100cc Sample	1.5 hours	455300005339/Gardco Gel Timer
Hardness	50 Shore A	455300006287/ASTM D2240
Glass Transition Temperature/Tg	-26 °C	453560822409 by DSC
Water Absorption	0.12 %	24 hr immersion 457561824543/ASTM D570
Peak Exotherm	38 °C after 18 minutes for 50 mL sample	455300005593 by Type K thermocouple

RESINLAB L.L.C. MAKES NO EXPRESS OR IMPLIED WARRANTIES OR MERCHANTABILITY, FITNESS OR OTHERWISE with respect to its products. In addition, while the information contained herein is believed to be reliable, no warranty is expressed or implied regarding the accuracy of the data or the results to be obtained from the use thereof. All recommendations or suggestions for use are made without guarantee inasmuch as conditions of use are beyond our control. The properties given are typical values and are not intended for use in preparing specifications. Users should make their own purposes. Page 1 of 3



TECHNICAL DATA SHEET UR3010RC Clear

Revision date: 2/25/2022

N109 W13300 ELLSWORTH DRIVE GERMANTOWN, WI 53022 262-253-5900 FAX 262-253-5919

Property:	Value:	Test Method or Source:
Tensile Properties:		4535601224470/ASTM D638
Strength	260 psi	
Elongation	93 %	
Modulus	440 psi	
Flame Resistance	Passes with HB Rating @ 6.0 mm	45376013225560/UL94HB
Tested at ResinLab, not UL Certified		
Volume Resistivity	2.25 x 10 ¹⁴ ohm-cm	455300006612/ASTM D257
		@ 23 °C @ 52 %RH
Surface Resistivity	1.70 x 10 ¹⁴ ohm/sq	455300006612/ASTM D257
		@ 23 °C @ 52 %RH
Dielectric Constant & Dissipation Factor		455300006513/ASTM D150
@ 100 Hz	6.5, 0.103	
@ 100 kHz	3.6, 0.095	
Coefficient of Thermal Expansion by TM	Α	455300005340/ASTM E831 TMA, 5 °C/min
above Tg	239 ppm/°C	
Operating Temperature Range	-40 to 125 °C**	
Relative Thermal Index (RTI)	50 °C	UL746B, Table 7.1
		Generic Value Based on Composition

* Asterisk denotes values considered typical to associated resin systems or extrapolated from other test results.

** Operating Temperature Range is based on average design requirements and is not intended as a guarantee of suitability for all applications operating at that temperature.

*** This TDS contains values that have been updated. The values reported in this technical data sheet are typical values of the product, and are highly dependent on test conditions and methodology. We actively seek the most precise and accurate ways to measure and interpret performance of our products, and to update estimated values with measured values. The formula has not been revised or changed in any way. Although the values on paper have changed, you can expect the same performance of the product.

INSTRUCTIONS:

- 1. Bring to room temperature prior to use.
- Cartridge format: Mixer should be attached keeping the cartridge vertical and any air pocket purged this way. After the mixer contains material, the mixer tip can be dropped to dispense pre-bleed amount. Attach a new static mixer with each cartridge, then pre-bleed the first 3 inches of dispensed material or until a uniform color is obtained. Maintain adequate velocity during dispensing to ensure complete mixing.
- 3. Bulk format: stir until homogeneous weigh and mix parts A and B accurately and thoroughly, scraping sides of container often. Do not pour from mixing container, transfer to a new container as residual unmixed material may cause a tacky spot on the surface of the casting. Maintain adequate velocity during dispensing to ensure complete mixing.
- 4. Allow to cure undisturbed until product is fully gelled or tack-free to the touch.

RESINLAB L.L.C. MAKES NO EXPRESS OR IMPLIED WARRANTIES OR MERCHANTABILITY, FITNESS OR OTHERWISE with respect to its products. In addition, while the information contained herein is believed to be reliable, no warranty is expressed or implied regarding the accuracy of the data or the results to be obtained from the use thereof. All recommendations or suggestions for use are made without guarantee inasmuch as conditions of use are beyond our control. The properties given are typical values and are not intended for use in preparing specifications. Users should make their own test to determine the suitability of this product for their own purposes. Page 2 of 3



Revision date: 2/25/2022

N109 W13300 ELLSWORTH DRIVE GERMANTOWN, WI 53022 262-253-5900 FAX 262-253-5919

5. Clean up uncured resin with suitable organic solvent such as MEK or acetone.

SHELF LIFE AND STORAGE:

6 months at 25 °C Bulk. 12 months at 25 °C in cartridges that are foil bagged and desiccant packed. Specialty packaging may be less.

Isocyanates are sensitive to moisture and should be kept in their original container or in a volume tank under dry nitrogen blanketing. Many isocyanates are prone to dimerization, the formation of a white precipitate. Products with minor amounts of this precipitate normally cure to full properties. Storage at 20 - 30 °C (68 °F to 86 °F) is recommended to ensure full shelf life.