

TECHNICAL DATA SHEET

EP750LV Clear

Revision date: 4/3/2023

N109 W13300 ELLSWORTH DRIVE GERMANTOWN, WI 53022 262-253-5900 FAX 262-253-5919

DESCRIPTION:

Resinlab® EP750LV Clear is a food grade two-part unfilled epoxy structural adhesive designed for bonding applications. It is comprised of raw materials listed in the Code of Federal Regulations Title 21:175.105 (Adhesives) and 175.300 (Resinous and polymeric coatings)¹. EP750LV Clear can be used safely in applications where there is direct or indirect food contact. This includes coatings intended for use in producing, manufacturing, packing, processing, preparing, treating, packaging, transporting or holding food.

EP750LV Clear cures completely at room temperature to a tough, semi rigid polymer. It has good wetting and adhesion to most surfaces and has a free-flowing viscosity. It has very good resistance to water, acids and bases and most organic solvents.

EP750LV Clear was formulated to a 2A:1B volume mix ratio for use in side by side dispensing cartridges and meter/mix and dispense equipment. *EP750LV Clear* will reach full cure at room temperature within 24-72 hours. Cure time can be accelerated by the application of heat. Times and temperatures from 2 hours at 65 °C to 15 minutes at 100 °C are typical for most applications. Time to heat substrate must be taken into account. Cooler temperatures will extend work time and increase cure times.

EP750LV Clear is suitable for use in medical device assembly, it has been tested and is proven non-toxic per ISO10993-5. Manufacturers should test their own finished product for biocompatibility. Certificates of compliance are available upon request.

TYPICAL PROPERTIES:

All properties given are at 25 °C unless otherwise noted.

Property:	Value:	Test Method or Source:
Color	Clear	Visual
Mix Ratio	Part A to Part B	Calculated
Mix Ratio by weight	2.44 to 1	
Mix Ratio by volume	2 to 1	
Cure Schedule	24-72 hrs @ 25 °C	
	2 hrs @ 65 °C	
	15 min @ 100 °C	
	10 min @ 110 °C	
Viscosity - Part A	14,500 cP	TA HR20 Rheometer 25mm parallel plate @
Viscosity - Part B	400 cP	1/s DCV6100723
Viscosity - Mixed	2,500 cP	
Specific Gravity - Part A	1.16	Calculated
Specific Gravity - Part B	0.95	
Specific Gravity - Mixed	1.09	



TECHNICAL DATA SHEET

EP750LV Clear

Revision date: 4/3/2023

N109 W13300 ELLSWORTH DRIVE GERMANTOWN, WI 53022 262-253-5900 FAX 262-253-5919

Property:	Value:	Test Method or Source:
Pot Life defined as the time it takes for	45 minutes	TA HR20 Rheometer parallel plate 25mm @
initial mixed viscosity to double		1/s DCV6100723
Gel Time 100cc Sample	4 hours	455300005339/Gardco Hot Pot Gel Timer
Hardness	80 Shore D	455300006287/ASTM D2240
Glass Transition Temperature/Tg	74 °C	453560822409 by DSC
Water Absorption	0.18 %	24 hr immersion 457561824543/ASTM D570
Tensile Properties:		4535601224470/ASTM D638
Strength	7,000 psi	
Elongation	3 – 5 %	
Modulus	315,000 psi	
Lap Shear Strength		4535601224468/ASTM D1002
0.010" Bond Line, Al to Al	1,500 psi	
Compressive Properties:		4535601224467/ASTM D695
Yield Strength	9,000 psi	
Ultimate Strength	11,000 psi	
Modulus	450,000 psi	
Electrical Resistivity:		455300006612/ASTM D257
Volume	8 x 10 ¹⁴ ohm-cm (@ 25 °C)*	
Surface	1 x 10 ¹⁵ ohm/sq (@ 50%RH)*	
Dielectric Constant & Dissipation Factor:		455300006513/ASTM D150
@ 100 Hz	3.2, 0.01*	
@ 100 kHz	3.1, 0.01*	
AC Dielectric Strength	17 kV/mm *	DCV6101609; ASTM D149 Method A,
		immersed in ASTM D3487 Type II Oil
Coefficient of Thermal Expansion by TMA:		455300005340/ASTM E831 TMA, 5 °C/min
below Tg	75 ppm/°C	
above Tg	203 ppm/°C	
Biocompatibility Biological Evaluation of	Passes ISO 10993-5	MEM Elution Test
Medical Devices		
Operating Temperature Range	-55 to 150 °C**	
Relative Thermal Index (RTI)	90 °C	UL746B, Table 7.1
		Generic Value Based on Composition

^{*} Asterisk denotes values considered typical to associated resin systems or extrapolated from other test results.

^{**} Operating Temperature Range is based on average design requirements and is not intended as a guarantee of suitability for all applications operating at that temperature.

^{***} This TDS contains values that have been updated. The values reported in this technical data sheet are typical values of the product, and are highly dependent on test conditions and methodology. We actively seek the most precise and accurate ways to measure and interpret performance of our products, and to update estimated values with measured values. The formula has not been revised or changed in any way. Although the values on paper have changed, you can expect the same performance of the product.



TECHNICAL DATA SHEET

EP750LV Clear

Revision date: 4/3/2023

N109 W13300 ELLSWORTH DRIVE GERMANTOWN, WI 53022 262-253-5900 FAX 262-253-5919

INSTRUCTIONS:

- 1. Bring to room temperature prior to use.
- 2. Cartridge format: Mixer should be attached keeping the cartridge vertical and any air pocket purged this way. After the mixer contains material, the mixer tip can be dropped to dispense pre-bleed amount. Attach a new static mixer with each cartridge, then pre-bleed the first 3 inches of dispensed material or until a uniform color is obtained. Maintain adequate velocity during dispensing to ensure complete mixing.
- 3. Bulk format: stir until homogeneous weigh and mix parts A and B accurately and thoroughly, scraping sides of container often. Do not pour from mixing container, transfer to a new container as residual unmixed material may cause a tacky spot on the surface of the casting. Maintain adequate velocity during dispensing to ensure complete mixing.
- 4. Clean up uncured resin with suitable organic solvent such as MEK or acetone.
- 5. Allow to cure undisturbed until product is fully gelled or tack-free to the touch.

SHELF LIFE AND STORAGE:

12 months at 25 °C. Specialty packaging may be less.

Many epoxy resin systems are prone to crystallization as epoxy resin is a super-cooled fluid. This condition may give the product a gritty or grainy appearance (or hazy in clear products). Products in this state will not usually cure to normal and expected properties. In extreme cases it may appear solid and cured. Fluctuating temperatures (within 5 to 50 °C) aggravate this phenomenon. Heating the individual component to 50 to 60 °C while stirring can usually restore products to original state. Storage at 25 +/- 10 °C is optimum for most products.

¹The FDA does not approve products for direct or indirect food contact. This formula is comprised solely of materials listed under these regulations and can be "used in accordance with" or "used in compliance with CFR Title 21 175.300". It the responsibility of the end user to determine if this conformance is appropriate to the final product's end uses. The specific requirements of an application for food, drug or potable water are detailed by the FDA, USDA (food plants) or NSF (water).